

Anti-Gemin 2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O14893

Item Number: CM0002888



Introduction

Rabbit monoclonal antibody targeting human Gemin 2 (SIP1), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for studying snRNP assembly and pre-mRNA splicing.

[Learn More](#)

Parameter	Value
Product Name	Gemin 2 Rabbit mAb
Remarks/Alias	SIP1; SIP1-delta; Gemin 2
Species	Human
Gene ID (Human)	8487
Gene ID	8487
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Gemin 2 (O14893).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14893
Protein Weight	32kDa
Shipping	Ice bag