


Anti-Gemin8 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nwz8

Item Number: CM0013291



Introduction

GEMIN8 rabbit polyclonal antibody validated for WB and ELISA. Designed for human GEMIN8, a gem-associated protein involved in snRNP assembly and pre-mRNA splicing. - Q9NWZ8

[Learn More](#)

Parameter	Value
Product Name	GEMIN8 Rabbit pAb
Remarks/Alias	FAM51A1; GEMIN8
Species	Human
GeneID	54960
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-242 of human GEMIN8 (NP_060326.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NWZ8
Protein Weight	29kDa
Shipping	Ice bag