

Anti-Gen1 Rabbit Polyclonal Antibody For Wb And Elisa - Q17Rs7

Item Number: CM0016856



Introduction

Rabbit polyclonal antibody against human GEN1, validated for WB and ELISA. Targets flap endonuclease GEN homolog 1 (Q17RS7), involved in Holliday junction resolution and DNA repair. Suitable for human and mouse samples.

[Learn More](#)

Parameter	Value
Product Name	GEN1 Rabbit pAb
Remarks/Alias	Gen; GEN1
Species	Human
Gene ID	348654
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q17RS7
Protein Weight	103kDa
Shipping	Ice bag