


Anti-Gephyrin Monoclonal Antibody For Wb, Elisa - Q9Nqx3

Item Number: CM0005561



Introduction

Anti-Gephyrin Monoclonal Antibody validated for WB, ELISA in Human, Mouse, Rat. Targets inhibitory synapse organizer and molybdenum cofactor biosynthesis enzyme. For neuroscience and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Gephyrin Rabbit mAb
Remarks/Alias	GPH; GEPH; HKPX1; GPHRYN; MOCODC; Gephyrin
Species	Human
GeneID (Human)	10243
GeneID	10243
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQX3
Protein Weight	80kDa/83kDa
Shipping	Ice bag