


## Anti-Gephyrin Polyclonal Antibody For Wb And Elisa - Q9Nqx3

Item Number: CM0026742



### Introduction

Rabbit polyclonal antibody against human Gephyrin (GPHN) for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for investigating inhibitory synapse anchoring and molybdenum cofactor biosynthesis.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | Gephyrin Rabbit pAb                                                                                                   |
| Remarks/Alias    | GPH; GEPH; HKPX1; GPHRYN; MOCODC; Gephyrin                                                                            |
| Species          | Human                                                                                                                 |
| GeneID (Human)   | 10243                                                                                                                 |
| GeneID           | 10243                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 520-769 of human Gephyrin (NP_065857.1) |
| Source           | Rabbit                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                 |
| Application      | WB, ELISA                                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                     |
| SWISS            | Q9NQX3                                                                                                                |
| Protein Weight   | 80kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |