


Anti-Get4 Rabbit Polyclonal Antibody For Wb, Elisa - Q7L5D6

Item Number: CM0020906



Introduction

High-quality rabbit polyclonal antibody against Golgi to ER traffic protein 4 homolog (GET4), validated for WB and ELISA. Reacts with human, mouse, and rat. Ideal for protein quality control and tail-anchored protein targeting research.

[Learn More](#)

Parameter	Value
Product Name	GET4 Rabbit pAb
Remarks/Alias	CEE; CDG2Y; TRC35; CGI-20; C7orf20; GET4
Species	Human
GeneID (human)	51608
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-327 of human GET4 (NP_057033.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7L5D6
Protein Weight	37kDa
Shipping	Ice bag