

Anti-Gf14-C Rabbit Polyclonal Antibody For Wb And Elisa - Q6Zkc0

Item Number: CM0022294



Introduction

Rabbit polyclonal antibody against rice GF14-C (14-3-3-like protein G-box factor 14-3-3 homolog C). Validated for Western blot and ELISA applications. Suitable for plant research and diagnostics. Swiss-Prot Q6ZKC0.

[Learn More](#)

Parameter	Value
Product Name	GF14-C Rabbit pAb
Remarks/Alias	GF14C; OJ1124_B05.7; GF14-C
Species	Oryza sativa
Gene ID (Oryza sativa)	4345634
Immunogen	Recombinant protein of Oryza sativa GF14-C
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q6ZKC0
Protein Weight	71kDa
Shipping	Ice bag