

Anti-Gfap Monoclonal Antibody (Conjugate) For Immunofluorescence - P14136

Item Number: CM0009110



Introduction

Conjugated GFAP Rabbit Monoclonal Antibody validated for immunofluorescence (IF-P). Reacts with human, mouse, and rat GFAP, a 50 kDa astrocyte marker. Suitable for neuroscience research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Conjugated GFAP Rabbit mAb
Remarks/Alias	ALXDRD
Species	Human
GeneID (human)	2670
GeneID	2670
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GFAP (P14136).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P
Cross Reactivity	Human, Mouse, Rat
SWISS	P14136
Protein Weight	50kDa
Shipping	Ice bag