

Anti-Gfap Monoclonal Antibody For Wb, If-F, If-P, Mihc, Elisa - P14136

Item Number: CM0003202



Introduction

Anti-GFAP rabbit monoclonal antibody validated for WB, IF-F, IF-P, mIHC, ELISA. Cross-reacts with human, mouse, rat GFAP (~50 kDa). Ideal for neuroscience and astrocyte research.

[Learn More](#)

Parameter	Value
Product Name	GFAP Rabbit mAb
Remarks/Alias	ALXDRD; GFAP
Species	Human
GeneID (Human)	2670
GeneID	2670
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, IF-P, mIHC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14136
Protein Weight	50 kDa
Shipping	Ice bag