

Anti-Gfap Monoclonal Antibody For Wb, If, Ihc-P, Elisa - P14136

Item Number: CM0008822



Introduction

GFAP Rabbit mAb (P14136) is a rabbit monoclonal antibody validated for WB, IF/ICC, IF-P, IHC-P, ELISA. It detects human GFAP and cross-reacts with mouse and rat. ~50 kDa.

[Learn More](#)

Parameter	Value
Product Name	GFAP Rabbit mAb
Remarks/Alias	ALXDRD
Species	Human
GeneID (Human)	2670
GeneID	2670
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-383 of human GFAP (NP_002046.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P14136
Protein Weight	50kDa
Shipping	Ice bag