

Anti-Gfap Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P14136

Item Number: CM0004981



Introduction

Mouse monoclonal GFAP antibody validated for WB, IF-P, IHC-P, and ELISA in human, mouse, and rat. Essential for astrocyte detection in CNS research, with one-stop sourcing support from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	GFAP Mouse mAb
Remarks/Alias	ALXDRD
Species	Human
GeneID (Human)	2670
GeneID	2670
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 60-383 of human GFAP (NP_002046.1).
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14136
Protein Weight	50 kDa
Shipping	Ice bag