

Anti-Gfap Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Ihc-P, And Elisa - P14136

Item Number: CM0017567



Introduction

High-quality rabbit polyclonal antibody against human GFAP, validated for WB, IF/ICC, IF-P, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for astrocyte marker studies in neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	GFAP Rabbit pAb
Remarks/Alias	ALXDRD; GFAP
Species	Human
GeneID (Human)	2670
GeneID	2670
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-75 of human GFAP (NP_002046.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14136
Protein Weight	50kDa
Shipping	Ice bag