


Anti-Gfer Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P55789

Item Number: CM0024752



Introduction

High-quality rabbit polyclonal antibody against human GFER (ALR/ERV1), validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat GFER. Ideal for mitochondrial redox biology and liver regeneration research.

[Learn More](#)

Parameter	Value
Product Name	GFER Rabbit pAb
Remarks/Alias	ALR; HPO; HSS; ERV1; HPO1; HPO2; HERV1; MMCHD; MPMCD; GFER
Species	Human
GeneID (Human)	2671
GeneID	2671
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 81-205 of human GFER (NP_005253.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55789
Protein Weight	23kDa
Shipping	Ice bag