


Anti-Gfm1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Rp9

Item Number: CM0026985



Introduction

Rabbit polyclonal antibody targeting human mitochondrial elongation factor G (GFM1). Validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for mitochondrial translation and disease research.

[Learn More](#)

Parameter	Value
Product Name	GFM1 Rabbit pAb
Remarks/Alias	EFG; GFM; EFG1; EFGM; EGF1; hEFG1; COXPD1; mtEF-G1; GFM1
Species	Human
GeneID (Human)	85476
GeneID	85476
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 482-751 of human GFM1 (NP_079272.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96RP9
Protein Weight	83kDa
Shipping	Ice bag