

Anti-Gfpt1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q06210

Item Number: CM0002584



Introduction

GFPT1 Rabbit mAb validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat GFPT1 (~79 kDa). Essential for hexosamine pathway and glycosylation studies. Detects a key enzyme in UDP-GlcNAc synthesis.

[Learn More](#)

Parameter	Value
Product Name	GFPT1 Rabbit mAb
Remarks/Alias	GFA; GFAT; GFPT; MSLG; CMS12; GFAT1; CMSTA1; GFAT 1; GFAT1m; GFPT1L; GFPT1
Species	Human
Gene ID (Human)	2673
Gene ID	2673
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 600-699 of human GFPT1 (Q06210)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06210
Protein Weight	79 kDa
Shipping	Ice bag