

Anti-Gfpt1 Rabbit Polyclonal Antibody (Wb, Elisa) - Q06210

Item Number: CM0025375



Introduction

High-quality rabbit polyclonal antibody against GFPT1, validated for Western blotting (WB) and ELISA. Suitable for human, mouse, and rat samples. Target hexosamine pathway enzyme with a molecular weight of 79kDa.

[Learn More](#)

Parameter	Value
Product Name	GFPT1 Rabbit pAb
Remarks/Alias	GFA; GFAT; GFPT; MSLG; CMS12; GFAT1; CMSTA1; GFAT 1; GFAT1m; GFPT1L; GFPT1
Species	Human
Gene ID (Human)	2673
Gene ID	2673
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human GFPT1 (NP_002047.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06210
Protein Weight	79kDa
Shipping	Ice bag