

Anti-Gfpt2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O94808

Item Number: CM0024022

Introduction

High-quality rabbit polyclonal antibody targeting human GFPT2 (SWISS: O94808). Validated for WB, IF/ICC, and ELISA; cross-reacts with mouse and rat. Suitable for hexosamine biosynthesis pathway and protein glycosylation research.

[Learn More](#)

Parameter	Value
Product Name	GFPT2 Rabbit pAb
Remarks/Alias	GFAT; GFAT2; GFAT 2; GFPT2
Species	Human
Gene ID	9945
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-260 of human GFPT2 (NP_005101.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94808
Protein Weight	77kDa
Shipping	Ice bag