


Anti-Gfra1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P56159

Item Number: CM0029101



Introduction

GFRA1 rabbit polyclonal antibody validated for Western blot, IHC-P, IF/ICC, and ELISA. Detects GFRA1 in human, mouse, and rat samples. Ideal for neurotrophic signaling and dopaminergic neuron research.

[Learn More](#)

Parameter	Specification
Product Name	GFRA1 Rabbit pAb
Remarks/Alias	GNFR; RET1L; RETL1; RHDA4; TRNR1; GDNFRA; GFRalpha-1; GFR-ALPHA-1; GDNFR-alpha-1; GFRA1
Species	Human
GeneID	2674
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 345-420 of human GFRA1 (NP_005255.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56159
Protein Weight	51kDa
Shipping	Ice bag