


Anti-Gfra2 Rabbit Polyclonal Antibody For Wb And Elisa - O00451

Item Number: CM0031030



Introduction

Affinity-purified rabbit polyclonal antibody against human GFRA2, validated for Western Blot and ELISA in human, mouse, and rat. Targets neurturin receptor alpha for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	GFRA2 Rabbit pAb
Remarks/Alias	NTNRA; RETL2; TRNR2; GDNFRB; NRTNR-ALPHA; GFRA2
Species	Human
GeneID (Human)	2675
GeneID	2675
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 360-464 of human GFRA2 (NP_001486.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00451
Protein Weight	52kDa
Shipping	Ice bag