

Anti-Gfra3 Rabbit Polyclonal Antibody For Wb And Elisa - O60609

Item Number: CM0030039



Introduction

Affinity-purified rabbit polyclonal antibody targeting human GFRA3 (GDNF family receptor alpha-3). Suitable for WB and ELISA; cross-reacts with human and mouse. Ideal for neurotrophic factor signaling and peripheral neuron research.

[Learn More](#)

Parameter	Value
Product Name	GFRA3 Rabbit pAb
Remarks/Alias	GDNFR3; GFRA3
Species	Human
GeneID (Human)	2676
GeneID	2676
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 32-240 of human GFRA3 (NP_001487.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60609
Protein Weight	45kDa
Shipping	Ice bag