


Anti-Gga2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ujy4

Item Number: CM0014158



Introduction

GGA2 Rabbit polyclonal antibody validated for Western blot, immunofluorescence/immunocytochemistry, and ELISA. Detects endogenous GGA2 in human and mouse samples. Ideal for studying protein trafficking and Golgi-endosome dynamics.

[Learn More](#)

Parameter	Value
Product Name	GGA2 Rabbit pAb
Remarks/Alias	VEAR; GGA2
Species	Human
GeneID (Human)	23062
GeneID	23062
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human GGA2 (NP_055859.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UJY4
Protein Weight	67kDa
Shipping	Ice bag