

Anti-Ggct Polyclonal Antibody For Wb, Elisa - 075223

Item Number: CM0031017



Introduction

Rabbit polyclonal antibody targeting human gamma-glutamylcyclotransferase (GGCT), validated for WB and ELISA. Ideal for studies on glutathione homeostasis and apoptosis. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	GGCT Rabbit pAb
Remarks/Alias	GGC; GCTG; CRF21; C7orf24; GGCT
Species	Human
GeneID (Human)	79017
GeneID	79017
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-188 of human GGCT (NP_076956.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	075223
Protein Weight	21kDa
Shipping	Ice bag