

Anti-Ggcx Rabbit Polyclonal Antibody For Wb And Elisa - P38435

Item Number: CM0021807



Introduction

GGCX Rabbit Polyclonal Antibody (P38435) validated for Western blot (WB) and ELISA. Detects human and mouse vitamin K-dependent gamma-carboxylase, essential for coagulation and bone metabolism research.

[Learn More](#)

Parameter	Value
Product Name	GGCX Rabbit pAb
Remarks/Alias	VKCFD1; GGcX
Species	Human
Gene ID (Human)	2677
Gene ID	2677
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 479-758 of human GGcX (NP_000812.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P38435
Protein Weight	88kDa
Shipping	Ice bag