


Anti-Ggh Polyclonal Antibody For Wb, Elisa - Q92820

Item Number: CM0030486



Introduction

Rabbit anti-human GGH (gamma-glutamyl hydrolase) polyclonal antibody validated for WB and ELISA with cross-reactivity in human and mouse. Ideal for folate metabolism and antifolate research. UniProt Q92820.

[Learn More](#)

Parameter	Value
Product Name	GGH Rabbit pAb
Remarks/Alias	GH; GATD10; GGH
Species	Human
GeneID (human)	8836
GeneID	8836
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-230 of human GGH (NP_003869.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92820
Protein Weight	36kDa
Shipping	Ice bag