


Anti-Ggps1 Polyclonal Antibody For Wb, Elisa - 095749

Item Number: CM0019186



Introduction

Rabbit polyclonal antibody targeting GGPS1, suitable for WB and ELISA; cross-reacts with human, mouse, rat. Ideal for research on geranylgeranyl pyrophosphate synthesis and related diseases.

[Learn More](#)

Parameter	Value
Product Name	GGPS1 Rabbit pAb
Remarks/Alias	GGPPS; MDHLO; GGPPS1; MUDHLOV; GGPS1
Species	Human
GeneID (Human)	9453
GeneID	9453
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 164-300 of human GGPS1 (NP_001032354.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95749
Protein Weight	35kDa
Shipping	Ice bag