

## Anti-Ggt5 Rabbit Polyclonal Antibody For Wb And Elisa - P36269

Item Number: CM0016593

### Introduction

Validated for WB and ELISA, this rabbit polyclonal antibody detects human and mouse GGT5, with an immunogen corresponding to amino acids 388-587. Suitable for glutathione metabolism and inflammation studies.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | GGT5 Rabbit pAb                                                                                                       |
| Remarks/Alias    | GGL; GGT 5; GGTLA1; GGT-REL; GGT5                                                                                     |
| Species          | Human                                                                                                                 |
| GeneID (Human)   | 2687                                                                                                                  |
| GeneID           | 2687                                                                                                                  |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 388-587 of human GGT5 (NP_001093251.1). |
| Source           | Rabbit                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                 |
| Application      | WB, ELISA                                                                                                             |
| Cross Reactivity | Human, Mouse                                                                                                          |
| SWISS            | P36269                                                                                                                |
| Protein Weight   | 62kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |