


Anti- γ -Catenin Monoclonal Antibody For Wb, Elisa - P14923

Item Number: CM0006009



Introduction

Rabbit monoclonal antibody targeting γ -Catenin (JUP). Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for studying cell adhesion and junctional complexes.

[Learn More](#)

Parameter	Value
Product Name	γ -Catenin Rabbit mAb
Remarks/Alias	PG; DP3; PDGB; PKGB; CTNNG; DP113; γ -Catenin
Species	Human
GeneID (Human)	3728
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 15-115 of human γ -Catenin (P14923)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14923
Protein Weight	82kDa
Shipping	Ice bag