


## Anti- $\Gamma$ -Catenin Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P14923

Item Number: CM0013704



### Introduction

Rabbit polyclonal antibody recognizing human and mouse  $\gamma$ -catenin (junction plakoglobin). Validated for WB, IP, and ELISA, detects an 82 kDa protein. Ideal for cell adhesion, cardiomyopathy, and dermatology studies.

[Learn More](#)

| Parameter        | Value                                                                                                                        |
|------------------|------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | $\gamma$ -Catenin Rabbit pAb                                                                                                 |
| Remarks/Alias    | PG, DP3, PDGB, PKGB, CTNNG, DPIII, $\gamma$ -Catenin                                                                         |
| Species          | Human                                                                                                                        |
| GeneID (Human)   | 3728                                                                                                                         |
| GeneID           | 3728                                                                                                                         |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human $\gamma$ -Catenin (NP_002221.1) |
| Source           | Rabbit                                                                                                                       |
| Category         | Polyclonal Antibodies                                                                                                        |
| Application      | WB, IP, ELISA                                                                                                                |
| Cross Reactivity | Human, Mouse                                                                                                                 |
| SWISS            | P14923                                                                                                                       |
| Protein Weight   | 82kDa                                                                                                                        |
| Shipping         | Ice bag                                                                                                                      |