

Anti- γ -Tubulin Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P23258

Item Number: CM0007339



Introduction

γ -Tubulin Rabbit Monoclonal Antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Recognizes gamma-tubulin, a key centrosome protein involved in microtubule nucleation and cell division. Ideal for cell cycle and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	γ -Tubulin Rabbit mAb
Remarks/Alias	TUBG; GCP-1; CDCBM4; TUBGCP1; γ -Tubulin
Species	Human
Gene ID (Human)	7283
Immunogen	Recombinant protein of human γ -Tubulin
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23258
Protein Weight	51kDa
Shipping	Ice bag