

Anti- γ -Tubulin Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P23258

Item Number: CM0006268



Introduction

Rabbit monoclonal antibody targeting γ -Tubulin (TUBG1), validated for WB, IHC-P, IF/ICC, ELISA. Reacts with Human, Mouse, Rat. Essential for microtubule nucleation and centrosome research.

[Learn More](#)

Parameter	Value
Product Name	γ -Tubulin Rabbit mAb
Remarks/Alias	TUBG; GCP-1; CDCBM4; TUBGCP1; γ -Tubulin
Species	Human
Gene ID (Human)	7283
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 231-451 of human γ -Tubulin (P23258)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23258
Protein Weight	51kDa
Shipping	Ice bag