

Anti- γ -Tubulin Rabbit Polyclonal Antibody For Wb, Elisa - P23258

Item Number: CM0029482



Introduction

γ -Tubulin Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with Human, Mouse, Rat. Targets gamma-tubulin, a key component of the γ -Tubulin Ring Complex (γ TuRC) involved in microtubule nucleation and centrosome function.

[Learn More](#)

Parameter	Specification
Product Name	γ -Tubulin Rabbit pAb
Remarks/Alias	TUBG; GCP-1; CDCBM4; TUBGCP1; γ -Tubulin
Species	Human
Gene ID	7283 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 382-451 of human γ -Tubulin (NP_001061.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23258
Protein Weight	51kDa
Shipping	Ice bag