

Anti-Ghd 7.1 Rabbit Polyclonal Antibody For Wb, Elisa - Q0D3B6

Item Number: CM0021514

Introduction

Rabbit polyclonal antibody targeting *Oryza sativa* Ghd 7.1 (PRR37), a flowering time repressor. Validated for WB and ELISA, with cross-reactivity to rice. Useful for photoperiodic flowering and circadian clock studies.

[Learn More](#)

Parameter	Value
Product Name	Ghd 7.1 Rabbit pAb
Remarks/Alias	HD2; DTH7; PRR37; OsPRR37; Ghd 7.1
Species	<i>Oryza sativa</i>
GeneID	4344399
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of <i>Oryza sativa</i> Ghd 7.1 (Q0D3B6)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Oryza sativa</i>
SWISS	Q0D3B6
Protein Weight	80kDa
Shipping	Ice bag