

## Anti-Ghd7 Polyclonal Antibody For Wb, Elisa - E5Rqa1

Item Number: CM0023505

### Introduction

Rabbit polyclonal antibody against rice transcription factor GHD7, validated for WB and ELISA. Targets a key floral repressor under long-day conditions, suitable for flowering time regulation studies in *Oryza sativa*.

[Learn More](#)

| Parameter        | Value                                    |
|------------------|------------------------------------------|
| Product Name     | GHD7 Rabbit pAb                          |
| Remarks/Alias    | Ghd7; GHD7                               |
| Species          | Oryza sativa                             |
| Gene ID          | 107276161                                |
| Immunogen        | Recombinant protein of Oryza sativa GHD7 |
| Source           | Rabbit                                   |
| Category         | Polyclonal Antibodies                    |
| Application      | WB, ELISA                                |
| Cross Reactivity | Oryza sativa                             |
| SWISS            | E5RQA1                                   |
| Protein Weight   | 27kDa                                    |
| Shipping         | Ice bag                                  |