


Anti-Ghrh Monoclonal Antibody For Wb, Elisa - P01286

Item Number: CM0005710



Introduction

Rabbit monoclonal antibody targeting human GHRH (Somatoliberin). Validated for WB and ELISA, this antibody reacts with a 12kDa protein. Suitable for growth hormone secretion research.

[Learn More](#)

Parameter	Value
Product Name	GHRH Rabbit mAb
Remarks/Alias	GRF; INN; GHRF; GHRH
Species	Human
Gene ID (Human)	2691
Gene ID	2691
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human GHRH. (NP_066567.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01286
Protein Weight	12kDa
Shipping	Ice bag