


Anti-Ghrh Polyclonal Antibody For Wb, If/Icc, Elisa - P01286

Item Number: CM0028920



Introduction

Rabbit polyclonal anti-GHRH antibody validated for Western blot, IF/ICC, and ELISA. Reacts with human, mouse, rat GHRH (somatoliberin). Ideal for growth hormone secretion research.

[Learn More](#)

Parameter	Value
Product Name	GHRH Rabbit pAb
Remarks/Alias	GRF; INN; GHRF; GHRH
Species	Human
Gene ID (Human)	2691
Gene ID	2691
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human GHRH (NP_066567.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01286
Protein Weight	12kDa
Shipping	Ice bag