

Anti-Ghrhr Polyclonal Antibody For Wb, If-P, Elisa - Q02643

Item Number: CM0030643



Introduction

Rabbit polyclonal antibody targeting human GHRHR (growth hormone-releasing hormone receptor). Suitable for WB, IF-P, and ELISA. Cross-reacts with human, mouse, and rat samples. Supports research in growth hormone signaling and deficiency.

[Learn More](#)

Parameter	Value
Product Name	GHRHR Rabbit pAb
Remarks/Alias	GRFR; GHRFR; IGHD4; IGHD1B; GHRHR
Species	Human
Gene ID (Human)	2692
Gene ID	2692
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human GHRHR (NP_000814.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02643
Protein Weight	47kDa
Shipping	Ice bag