

Anti-Ghrl Polyclonal Antibody For Wb And Elisa - Q9Ubu3

Item Number: CM0018591



Introduction

High-quality rabbit polyclonal antibody against human GHRL (ghrelin). Validated for WB and ELISA. Ideal for studying appetite regulation, growth hormone release, and gastric function.

[Learn More](#)

Parameter	Value
Product Name	GHRL Rabbit pAb
Remarks/Alias	MTLRP; GHRL
Species	Human
GeneID (Human)	51738
GeneID	51738
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-117 of human GHRL (NP_057446.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UBU3
Protein Weight	13kDa
Shipping	Ice bag