

## Anti-Ghsr Rabbit Polyclonal Antibody For Wb, Elisa - Q92847

Item Number: CM0023714



### Introduction

Rabbit polyclonal antibody against human GHSR (ghrelin receptor). Suitable for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for research on appetite regulation and growth hormone secretion.

[Learn More](#)

| Parameter        | Value                                                                                              |
|------------------|----------------------------------------------------------------------------------------------------|
| Product Name     | GHSR Rabbit pAb                                                                                    |
| Remarks/Alias    | GHDP; GHSR                                                                                         |
| Species          | Human                                                                                              |
| Gene ID (Human)  | 2693                                                                                               |
| Gene ID          | 2693                                                                                               |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 1-100 of human GHSR (NP_940799.1) |
| Source           | Rabbit                                                                                             |
| Category         | Polyclonal Antibodies                                                                              |
| Application      | WB, ELISA                                                                                          |
| Cross Reactivity | Mouse, Rat                                                                                         |
| SWISS            | Q92847                                                                                             |
| Protein Weight   | 41kDa                                                                                              |
| Shipping         | Ice bag                                                                                            |