

Anti-Giantin Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q14789

Item Number: CM0021295

Introduction

Rabbit polyclonal antibody against human Giantin (GOLGB1, Q14789). Validated for WB, IF/ICC, ELISA. Targets a ~376 kDa Golgi membrane protein involved in intercisternal cross-bridge formation.

[Learn More](#)

Parameter	Value
Product Name	Giantin Rabbit pAb
Remarks/Alias	GCP; GCP372; GOLIM1; Giantin
Species	Human
GeneID	2804
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1123-1302 of human Giantin (NP_004478.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q14789
Protein Weight	376kDa
Shipping	Ice bag