

Anti-Gid1A Rabbit Polyclonal Antibody For Wb - Q9Maa7

Item Number: CM0024230



Introduction

Rabbit polyclonal antibody against Arabidopsis thaliana gibberellin receptor GID1A (ATCXE10). Validated for Western blot (WB) and ELISA. Ideal for studying GA signaling, plant growth and development. Immunogen: synthetic peptide within aa 50-150.

[Learn More](#)

Parameter	Value
Product Name	GID1A Rabbit pAb
Remarks/Alias	ATGID1A, GA INSENSITIVE DWARF1A, T12H1.8, T12H1_8, GID1A
Species	Arabidopsis thaliana
GeneID	819674
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of Arabidopsis thaliana GID1A (NP_187163.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
SWISS	Q9MAA7
Protein Weight	39kDa
Shipping	Ice bag