


Anti-Gif Rabbit Polyclonal Antibody For Wb, Elisa - P27352

Item Number: CM0021733



Introduction

High-quality anti-GIF rabbit polyclonal antibody validated for WB and ELISA, cross-reactive with human, mouse, and rat. Ideal for studying cobalamin absorption and intrinsic factor-related research.

[Learn More](#)

Parameter	Value
Product Name	GIF Rabbit pAb
Remarks/Alias	GIF; IF; IFMH; INF; TCN3; gastric intrinsic factor
Species	Human
GeneID (human)	2694
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-417 of human GIF (NP_005133.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27352
Protein Weight	42kDa/45kDa
Shipping	Ice bag