

Anti-Gigyf2 Rabbit Polyclonal Antibody For Wb, Elisa - Q6Y7W6

Item Number: CM0015494



Introduction

Anti-GIGYF2 rabbit polyclonal antibody validated for WB, ELISA. Reacts with Human, Mouse, Rat. Suitable for studies on translational repression, mRNA decay, and Parkinson's disease.

[Learn More](#)

Parameter	Value
Product Name	GIGYF2 Rabbit pAb
Remarks/Alias	GYF2; PERQ2; PERQ3; PARK11; TNRC15; GIGYF2
Species	Human
GeneID (human)	26058
GeneID	26058
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-500 of human GIGYF2 (NP_001096617.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6Y7W6
Protein Weight	150kDa
Shipping	Ice bag