


Anti-Gimap1 Polyclonal Antibody For Wb, Elisa - Q8Wwp7

Item Number: CM0022758



Introduction

High-quality rabbit polyclonal antibody targeting human GIMAP1, validated for WB and ELISA. Cross-reacts with rat. Essential for studying lymphocyte survival and GTPase biology in spleen and lymph nodes.

[Learn More](#)

Parameter	Value
Product Name	GIMAP1 Rabbit pAb
Remarks/Alias	IAN2; IMAP1; HIMAP1; IMAP38; GIMAP1
Species	Human
Gene ID (Human)	170575
Gene ID	170575
Immunogen	Recombinant protein corresponding to a sequence within amino acids 140-280 of human GIMAP1 (NP_570115.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q8WWP7
Protein Weight	34kDa
Shipping	Ice bag