

Anti-Gimap5 Rabbit Polyclonal Antibody For Wb, Elisa - Q96F15

Item Number: CM0030938



Introduction

High-quality rabbit polyclonal antibody against human GIMAP5, validated for WB and ELISA. Detects endogenous protein in human, mouse, rat samples. Essential for studies on T lymphocyte development and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	GIMAP5 Rabbit pAb
Remarks/Alias	IAN4; IAN5; IROD; IAN-5; IMAP3; NCPH2; HIMAP3; IAN4L1; GIMAP5
Species	Human
GeneID (Human)	55340
GeneID	55340
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 190-270 of human GIMAP5 (NP_060854.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96F15
Protein Weight	35kDa
Shipping	Ice bag