

Anti-Gins2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y248

Item Number: CM0028601



Introduction

High-quality rabbit polyclonal antibody against human GINS2 (PSF2). Validated for WB and ELISA. Ideal for DNA replication research, particularly the CMG helicase complex. Cross-reacts with human. Predicted molecular weight 21 kDa.

[Learn More](#)

Parameter	Value
Product Name	GINS2 Rabbit pAb
Remarks/Alias	PSF2; Pfs2; HSPC037; GINS2
Species	Human
GeneID	51659
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-185 of human GINS2 (NP_057179.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y248
Protein Weight	21kDa
Shipping	Ice bag