


Anti-Gins4 (Sld5) Rabbit Polyclonal Antibody For Wb, Elisa - Q9Brt9

Item Number: CM0030356



Introduction

Rabbit polyclonal antibody against human GINS4 (SLD5) for Western blot and ELISA. Detects human, mouse, and rat GINS4, a DNA replication complex protein. Ideal for cell cycle and replication research.

[Learn More](#)

Parameter	Value
Product Name	GINS4 Rabbit pAb
Remarks/Alias	SLD5, GINS4
Species	Human
Gene ID	84296
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BRT9
Protein Weight	98kDa
Shipping	Ice bag