

Anti-Gip Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P09681

Item Number: CM0030091



Introduction

Rabbit polyclonal antibody against human GIP (Gastric inhibitory polypeptide/Glucose-dependent insulinotropic polypeptide). Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Suitable for diabetes and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	GIP Rabbit pAb
Remarks/Alias	GIP
Species	Human
Gene ID (Human)	2695
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-153 of human GIP (NP_004114.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09681
Protein Weight	17kDa
Shipping	Ice bag