

Anti-Gipc1 Rabbit Polyclonal Antibody For Wb, Elisa - O14908

Item Number: CM0014781

Introduction

Rabbit polyclonal antibody targeting human GIPC1 (O14908), validated in WB and ELISA, with cross-reactivity to mouse and rat. Ideal for G protein-linked signaling studies.

[Learn More](#)

Parameter	Value
Product Name	GIPC1 Rabbit pAb
Remarks/Alias	NIP; GIPC; IIP-1; OPDM2; TIP-2; SEMCAP; C19orf3; Hs.6454; GLUT1CBP; RGS19IP1; SYNECTIN; SYNECTIIN; GIPC1
Species	Human
GeneID (Human)	10755
GeneID	10755
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-333 of human GIPC1 (NP_005707.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14908
Protein Weight	36kDa
Shipping	Ice bag