

Anti-Gipc2 Polyclonal Antibody For Wb And Elisa - Q8Tf65

Item Number: CM0013113



Introduction

Rabbit polyclonal antibody targeting human GIPC2 (PDZ domain-containing protein GIPC2, SEMCAP2) for Western blotting and ELISA applications. Recognizes a ~34 kDa cytoplasmic protein, down-regulated in kidney and colorectal tumors.

[Learn More](#)

Parameter	Value
Product Name	GIPC2 Rabbit pAb
Remarks/Alias	SEMCAP2; SEMCAP-2; GIPC2
Species	Human
GeneID (Human)	54810
GeneID	54810
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-315 of human GIPC2 (NP_060125.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8TF65
Protein Weight	34kDa
Shipping	Ice bag