

Anti-Girk2 (Kcnj6) Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - P48051

Item Number: CM0006420

Introduction

Rabbit monoclonal antibody against GIRK2 (KCNJ6), a G protein-activated inward rectifier potassium channel. Validated for WB, IF/ICC, IF-P, ELISA; reacts with human and rat. Useful for studying insulin secretion and ion channel-related diseases.

[Learn More](#)

Parameter	Specification
Product Name	GIRK2(KCNJ6) Rabbit mAb
Remarks/Alias	BIR1; GIRK2; KATP2; KCNJ7; KPLBS; GIRK-2; KATP-2; KIR3.2; hiGIRK2; GIRK2(KCNJ6)
Species	Human
GeneID	3763
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P48051
Protein Weight	48kDa
Shipping	Ice bag